medi GmbH & Co. KG
Medicusstraße 1
95448 Bayreuth
Germany
T +49 921 912-0
F +49 921 912-780
ortho@medi.de
www.medi.de

medi Australia Pty Ltd
 83 Fennell Street
 North Parramatta NSW 2151
 Australia
 T +61-2 9890 8696
 F +61-2 9890 8439
 sales@mediaustralia.com.au
 www.mediaustralia.com.au

medi Austria GmbH
 Adamgasse 16/7
 6020 Innsbruck
 Austria
 T +43 512 57 95 15
 F +43 512 57 95 15 45
 vertrieb@medi-austria.at
 www.medi-austria.at

medi Belgium NV
 Staatsbaan 77/0099
 3945 Ham
 Belgium
 T: + 32-11 24 25 60
 F: +32-11 24 25 64
 info@medibelgium.be
 www.medibelgium.be

medi Canada Inc / médi Canada Inc
 104-1375 Lionel-Boulet,
 Varennes, Québec,
 QC Canada J3X 1P7
 T +1 450-583-3317 / +1 800 361 3153
 F +1 888-583-6827
 service@medicanada.ca
 www.medicanada.ca

MAXIS a.s.,
 medi group company
 Slezská 2127/13
 120 00 Prague 2
 Czech Republic
 T: +420 571 633 510
 F: +420 571 616 271
 info@maxis-medica.com
 www.maxis-medica.com

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medi Trading (Shanghai) Co. Ltd.
 Room 1703, T1, CES WEST BUND CENTER,
 No. 277, Long Lan Road
 Xuhui District
 200232 Shanghai
 People's Republic of China
 T: +86-21 50582319
 F: +86-21 50582319

Medi Turk Ortopedi Medikal İc ve Dis Tic. AS
 Mustafa Kemal Mah. 2156 Sokak No: 16 /4-5
 06510 Çankaya Ankara
 Turkey
 T: +90 312 435 20 26
 F: +90 312 434 17 67
 info@medi-turk.com
 www.medi-turk.com

UK Responsible Person
 medi UK Ltd.
 Plough Lane
 Hereford HR4 0EL
 Great Britain
 T +44-1432 37 35 00
 F +44-1432 37 35 10
 enquiries@mediuk.co.uk
 www.mediuk.co.uk

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Important notes

This medical device **MD** is made only for multiple uses on a single patient. If it is used for treating more than one patient, the manufacturer's product liability will become invalid. If you experience excessive pain or any discomfort whilst wearing, remove the product and immediately contact your doctor or medical supply store. Only wear the product on unbroken skin or over injured skin that has been medically treated, not directly on injured or damaged skin, and only after previous consultation with a medical professional.

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medi Epico ROMs

Intended purpose

medi Epico ROMs is an elbow brace for mobilisation with immobilisation of the proximal radioulnar articulation.

Indications

All indications that require early functional mobilisation of the elbow joint with limitation of the range of motion and immobilisation of the proximal radioulnar articulation, such as::

- Surgical fixation of medial and/or lateral collateral ligaments
- Conservative therapy of dislocated elbows
- Surgery on dislocated elbow fractures
- Medial or lateral collateral ligament reconstructive surgery
- Surgically treated bicep or tricep tears
- Treatment of radial head fracture (nonsurgical, postoperative)
- Treatment following insertion of a prothesis
- Distal humeral fractures (nonsurgical, postoperative)

Contraindications

- Concomitant neurological damage, e.g. ulnar nerve
- Compartment syndrome

Risks / Side Effects

If the aids are too tight, this can cause local pressure sores or constriction of blood vessels or nerves. Therefore you should consult the doctor treating you before you use it if the following circumstances apply to you:

- Condition or injury of the skin in the area of application, particularly with infectious indications (excessive heat,

swelling or reddening)

- Paraesthesia and circulatory disorders (e.g. as seen in diabetes, varicose veins)
- Lymph drainage disorders – equally, swelling of soft tissue adjacent to the area of application which is not related to a condition

When wearing snug-fitting aids, this can result in localised skin inflammation or irritation that is due to mechanical irritation of the skin (particularly in conjunction with sweating) or due to the composition of the material.

Intended patient groups

Healthcare professionals and patients, including persons who play a supporting role in care, are the intended users for this device, provided that they have received briefing on this from healthcare professionals.

Patient target group: Healthcare professionals should provide care to adults and children, applying the available information on the measurements/sizes and required functions/indications, in line with the information provided by the manufacturer, and acting under their own responsibility.

Setting the flexion and extension limits

Your medi Epico ROMs brace is pre-installed at an extension and flexion angle of 0°. The ChronoLock device makes sure that the degree setting of the brace joint is not changed unintentionally (e.g. by hitting an object). For setting of the flexion and/or extension limit, always turn the ChronoLock wheel until the control window is located below the arrow of the grey fixation button. Now, press the grey button to set the flexion limits first, followed by setting the extension limit. For safety reasons, turn the ChronoLock wheel into a neutral position out of the

flexion and extension setting range.

Flexion limit: 0°-120°

Extension limit: 0°-120°

Fitting instructions

- Open all the straps and put the brace on so that its joint lies directly over the lateral epicondyle. You can adjust the length of the arm rest with the two screws on the forearm splint.
- First close the two lower arm belts (No. 1+2), then the upper arm belts (No. 3+4).

The pads can be trimmed to size along the guide lines, using scissors, if necessary. The straps can be shortened. To achieve this, loosen the fastening ends and shorten the straps as required. Then fasten the strap ends again so that they are secure.

- Guide one end of the neck strap through the ring on the arm cuff and fasten the hook and loop.
- Pass the other end over your shoulder, guide it through the ring on the hand rest, and fasten the strap.

The frame of the orthosis is made of aluminium and can be adjusted as required. Avoid making alterations close to the joint area, however.

Important notes

Adjust the extension and flexion limits on the instructions of your doctor only. Particularly during the first few days after traumatic events, the belts should be slackened often so as to avoid severe swelling.

Care instructions

Hook and loop fasteners should be closed for washing. Soap residues, lotions and ointments can cause skin irritation and material wear.

The cushion for the hand support can be removed for washing. To remove the cushion, loosen the screw of the shoulder strap with the attached Allen key.

- Wash the paddings and straps by hand, preferably using medi clean washing agent.
- Do not bleach.
- Leave to dry naturally.
- Do not iron.
- Do not dry clean.

With normal use, the hinge mechanism should not require any maintenance, but it can be lubricated with PTFE spray.

Storage instructions

Keep the product in a dry place and do not expose to direct sunlight.

Material composition

Aluminum, PU foam, cotton

Liability

The manufacturer's liability will become void if the product is not used as intended. Please also refer to the corresponding safety information and instructions in this manual.

Disposal

The product can be disposed of in the domestic waste.

Your medi Team

Wishes you a quick return to full fitness.

In the event of any complaints regarding the product such as damage to the fabric or a fault in the fit, please report to your specialist medical retailer directly. Only serious incidents which could lead to a significant deterioration in health or to death are to be reported to the manufacturer or the relevant authorities in the EU member state. The criteria for serious incidents are defined in Section 2, No. 65 of the Regulation (EU) 2017/745 (MDR).

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